



OPTIKA®
M I C R O S C O P E S
I T A L Y

Digital Microscope with Camera & Tablet

- Innovative digital technology for microscope imaging
- 4 Plan-Achro objectives 4X, 10X, 40X, 100XOil/Water

Optika TB Series

- X-LED3 illumination with white LED
- CMOS camera, 2.0 USB
- Continuous auto exposure & compressed video capture
- Ergonomic design for comfortable long-term use
- Detachable tablet PC for analysis on the run with Optika Vision Lite software

ONLY
\$2,500

CAT#: B-290TB



Prices are valid until December 31, 2023. Price changes and/or early termination of this promotion may occur at any time. FREE SHIPPING.
Shipments may be delayed if we are experiencing a high volume of orders. Terms and conditions apply.