




CRYOGENIC VIALS



- Temps as low as -196°C
- Self-standing star foot
- Leakproof silicone O-ring
- Resealable bags of 50
- Gamma sterilized
- Free locking rack
- RNase, DNase, pyrogen, & endotoxin free
- Autoclavable

SPECIFICATIONS	
Material (tube/cap):	LN ₂ grade polypropylene
Material (o-ring/washer):	Silicone
Temp Range:	-196°C to +121°C*
Vial Diameter:	12.5mm
*LN ₂ gas phase only; not for liquid phase	



CRYOGENIC VIALS						
Cat. No	Description	Cap Thread	Height (with cap)	Vol.	Qty	List Price
V3801	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	40mm	1mL	500/cs	\$131.00
V3802	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	47.3mm	2mL	500/cs	\$138.00
V3803	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	66.6mm	3mL	500/cs	\$175.50
V3805	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	90.7mm	5mL	500/cs	\$198.00
V4801	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	40.2mm	1mL	500/cs	\$131.00
V4802	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	47.5mm	2mL	500/cs	\$138.00
V4803	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	60.5mm	3mL	500/cs	\$175.50
V4805	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	91mm	5mL	500/cs	\$198.00
V5809-A*	Color Cap inserts for Cryogenic Vials, assorted colors (white, blue, green, red, yellow), 500/pk	—	—	—	500/pk	\$18.00
V5920	Locking Rack for 50 "star foot" vials, 210 × 100 × 25mm	—	—	—	4/pk	\$49.00

ACCESSORIES			
Cat. No	Description	Qty	List Price
R1060	Freezer Box, slip on lid, polypropylene, 81 × 2 ml, 130 × 30 × 45 mm	5/pk	\$42.00
R2281	Freezer Box, hinged lid, polypropylene, 81 × 2 ml, 130 × 30 × 53 mm	5/pk	\$59.88
R2200	Freezer Box, hinged lid, polypropylene, 100 × 2 ml, 130 × 30 × 53 mm	5/pk	\$59.88
R2781-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 81 × 2mL, 133 × 133 × 53mm	5/pk	\$17.70
R2700-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 100 × 2mL, 133 × 133 × 53mm	5/pk	\$18.42
R3781	FlipTop™ Freezer Box, hinged lid, cardboard, 81 × 3-4mL, 133 × 133 × 75mm	50/cs	\$148.73
R3700	FlipTop™ Freezer Box, hinged lid, cardboard, 100 × 3-4mL, 133 × 133 × 75mm	50/cs	\$148.73

* Assorted color pack, also available in individual colors

